

# LIPOZERO G150

## Aesthetic Cavitation

### 52 PROGRAMS

**9 CaviMax150 – 10 MF-Plus - 5 Adipose Mobilization NorAd-526**

**2 Man Adiposity - 8 Beauty – 8 Skin – 10 Health**

Ultrasonic CAVITATION is the most efficient methodology to fight Cellulite and to eliminate excess fat

#### MECHANISM:

- When adipose tissue is submitted to deep low frequency ultrasonic energy, a succession of cyclical compression and decompression which produce cavities containing gases and vapours.
- Thanks to this process, the cavities continue to expand and compress till they explode, destroying in the process the fat cell membrane.
- The fat contained in the adipocytes is freed and eliminated through blood circulation without damaging the surrounding tissues.

#### THE GLOBUS PROFESSIONAL CAVITATION FREQUENCIES

The cavitation effect is efficiently produced only with low frequencies, lower than 600 KHz, because they are the only one able to penetrate in depth to destroy the adipose tissue. Lower the frequency, deeper is the action of the ultrasonic wave.

- **150 KHz** for adiposity and Cellulite less thick and less deep (Professional)
- **526 KHz** to associate the adipose destruction with the production of Noradrenalin (NorAd-526) which helps to burn fat thanks to the increase of the metabolism. Specially used for the adiposity of the men abdomen
- **3 MHz** for superficial beauty treatments and for skin treatments.

#### MF PLUS - MULTIFREQUENCY TREATMENTS: GLOBUS WORLD FIRST

Scientific researches have demonstrated that the application of Cavitation with Frequency Modulation (alternation of different frequencies) produce better results because the association of different frequencies allows to act at various depth and allows to personalize the treatment according to the type of cellulite.

#### THERAPEUTIC PROGRAMS

LIPOZERO G150 is provided with a wide range of therapeutic programs for the treatment of the principal pathologies at a **1 MHz** frequency.

**NEW STYLISH CARRYING CASE:** LIPOZERO G150 is supplied with a convenient and modern aluminium carrying case to transport it easily and in a stylish way.

**52 PROGRAMS - 30 free programs storage:**

- 9 CAVITATION CAVIMAX 150 KHz
- 10 MF-PLUS — MULTIFREQUENCY TREATMENTS
- 5 ADIPOSE MOBILIZATION NorAd-526
- 2 MAN ADIPOSITTY
- 8 BEAUTY
- 8 SKIN
- 10 THERAPY



#### TECHNICAL FEATURES

- Power: 3 W/cm<sup>2</sup>
- Output: 1 - Head: 1 Ø mm42 - IP 65
- Power Supply: main 220V 50/60 Hz and rechargeable batteries (in option)
- Dimensions: 17 x 22 x 6 cm
- Backlight Display
- Emission: continued and pulsed
- Emission Frequency: **3 MHz, 1 MHz, 526 KHz e 150 KHz**
- CE Device
- Standard: EN 60601-1 EN 60601-1-2

#### EQUIPMENT

- 1 carrying case
- 1 ultrasound Lipozero G150
- 1 Multi-frequency Head (**3 e 1 MHz, 526 e 150 KHz**)
- 1 contact gel bottle
- 1 Power Pack
- 1 operating manual



# GLOBUS

# LIPOZERO G150

## Program List

### CAVIMAX 150 KHz

Lipolysis  
Lipolysis Addominal  
Cellulite  
Cellulite Oedematous  
Cellulite Fibroses  
Culotte de cheval  
Cellulite Oedematous Legs  
Cellulite Fibroses Legs  
Cellulite with orange skin effects

### MF-PLUS

Lipolysis  
Cellulite  
Cellulite Oedematous  
Cellulite Fibroses  
Culotte de cheval  
Cellulite with orange skin effects  
Man Adiposity abdomen  
Man Adiposity Hip  
Man Adiposity Lower limbs  
Man Adiposity Upper limbs

### ADIPOSE MOBILISATION NorAd-526

Lipolysis-NorAd-526  
Cellulite- NorAd-526  
Adiposity - NorAd-526  
NorAd-526 1  
NorAd-526 2

### MAN ADIPOSITY

Man Adiposity abdomen  
Man Adiposity Hip

### BEAUTY

Massage  
Modelling  
Modelling Legs  
Drainage  
Drainage Legs  
Drainage arms  
Stretch marks  
Active principal application

### SKIN

Acne  
Face Lifting  
Wrinkle Treatment  
Scar  
Couperose  
Body spots  
Face spots,  
Restructuring skin

### HEALTH

Sciatica  
Lumbago  
Arthritis  
Carpal tunnel  
Neuralgia  
Tendinitis  
Epycondilitis  
Periarthritis  
Sprain  
Confusions



# GLOBUS