

# LIPOZERO G39

## Professional Cavitation

### 64 PROGRAMS

**4 UltraCav39 - 6 CaviMax150 - 10 MF-Plus - 5 Adipose Mobilization  
NorAd-526 - 4 Man Adiposity - 10 Beauty - 10 Skin - 15 Health**

Ultrasonic CAVITATION is the most efficient methodology to fight Cellulite and to eliminate excess fat

#### MECHANISM:

- When adipose tissue is submitted to deep low frequency ultrasonic energy, a succession of cyclical compression and decompression which produce cavities containing gases and vapours.
- Thanks to this process, the cavities continue to expand and compress till they explode, destroying in the process the fat cell membrane.
- The fat contained in the adipocytes is freed and eliminated through blood circulation without damaging the surrounding tissues.

#### THE GLOBUS PROFESSIONAL CAVITATION FREQUENCIES

The cavitation effect is efficiently produced only with low frequencies, lower than 600 KHz, because they are the only one able to penetrate in depth to destroy the adipose tissue. Lower the frequency, deeper is the action of the ultrasonic wave.

- **39 KHz** for adiposity and Cellulite particularly thick and deep that need an strong action (Top Professional)
- **150 KHz** for adiposity and Cellulite less thick and less deep (Professional)
- **526 KHz** to associate the adipose destruction with the production of Noradrenalin (NorAd-526) which helps to burn fat thanks to the increase of the metabolism. Specially used for the adiposity of the men abdomen
- **3 MHz** for superficial beauty treatments and for skin treatments.

#### MF PLUS - MULTIFREQUENCY TREATMENTS: GLOBUS WORLD FIRST

Scientific researches have demonstrated that the application of Cavitation with Frequency Modulation (alternation of different frequencies) produce better results because the association of different frequencies allows to act at various depth and allows to personalize the treatment according to the type of cellulite.

#### THERAPEUTIC PROGRAMS

LIPOZERO G150 is provided with a wide range of therapeutic programs for the treatment of the principal pathologies at a **1 MHz** frequency.

**NEW STYLISH CARRYING CASE:** LIPOZERO G39 is supplied with a convenient and modern aluminium carrying case to transport it easily and in a stylish way.

#### 64 PROGRAMS 30 free programs storage:

- 4 CAVITATION ULTRACAV 39 KHz
- 6 CAVITAZIONE CAVIMAX 150 KHz
- 10 MF PLUS — MULTIFREQUENCY TREATMENTS
- 5 ADIPOSE MOBILIZATION NorAd-526
- 4 MAN ADIPOSITO
- 10 BEAUTY
- 10 SKIN
- 15 THERAPY



#### TECHNICAL FEATURES

Power: 3 W/cm<sup>2</sup>  
 Output: 1 - Head: 2 - 1x Ø mm50 & 1 x Ø mm42 - IP65  
 Power Supply: main 220V 50/60 Hz or rechargeable batteries (in option)  
 Dimensions: 17 x 22 x 6 cm  
 Backlight Display  
 Emission: continued and pulsed  
 Emission Frequency: 3 & 1 MHz, 526, 150 & 39 KHz  
 CE Device  
 Standard: EN 60601-1 EN 60601-1-2

#### EQUIPMENT

- 1 carrying case
- 1 ultrasound Lipozero G150
- 1 Multi-frequency Head (3 e 1 MHz, 526 e 150 KHz )
- 1 Head 39KHz
- 1 contact gel bottle
- 1 Power Pack
- 1 operating manual



# GLOBUS

# LIPOZERO G39

## Program List

### ULTRA-CAV 39 KHz

Total Body 3-6 cm  
Total Body > 6 cm  
Thigh-gluteus 3-6 cm  
Thigh-gluteus > 6 cm

### CAVI-MAX 150 KHz

Lipolysis  
Lipolysis Addominal  
Cellulite  
Cellulite Oedematous  
Cellulite Fibroses  
Cellulite with orange skin effects

### MF-PLUS

Lipolysis  
Cellulite  
Cellulite Oedematous  
Cellulite Fibroses  
Culotte de cheval  
Cellulite with orange skin effects  
Man Adiposity abdomen  
Man Adiposity Hip  
Man Adiposity Lower limbs  
Man Adiposity Upper limbs

### ADIPOSE MOBILISATION

#### NorAd-526

Lipolysis-NorAd-526  
Cellulite- NorAd-526  
Adiposity - NorAd-526  
NorAd-526 1  
NorAd-526 2

### MAN ADIPOSITY

Man Adiposity abdomen  
Man Adiposity Hip  
Man Adiposity maschile Lower Limbs  
Man Adiposity maschile Lower Limbs

### BEAUTY

Toning  
Massage  
Modelling  
Modelling Legs  
Modelling abdomen  
Drainage  
Drainage Legs  
Drainage arms  
Stretch marks  
Active principal application

### SKIN

Acne  
Face Lifting  
Wrinkle Treatment  
Scar  
Couperose  
Body spots  
Face spots,  
Restructuring skin  
Skin smoothing  
Photoaging

### HEALTH

Sciatica  
Lumbago  
Arthritis  
Carpal tunnel  
Neuralgia  
Cruralgy  
Brachialgy  
Adhesion cicatrization  
Borsitis & Capsulitis  
Tendinitis  
Epycondilitis  
Periarthritis  
Sprain  
Contusions  
Ankle pathology



# GLOBUS