

LIPOZERO G500

Aesthetic Cavitation

21 PROGRAMS

5 Adipose Mobilization NorAd 526 – 4 Cavitation 330 KHz – 6 Beauty – 6 Skin

Ultrasonic CAVITATION is the most efficient methodology to fight Cellulite and to eliminate excess fat

MECHANISM:

- When adipose tissue is submitted to deep low frequency ultrasonic energy, a succession of cyclical compression and decompression which produce cavities containing gases and vapours.
- Thanks to this process, the cavities continue to expand and compress til they explode, destroying in the process the fat cell membrane.
- The fat contained in the adipocytes is freed and eliminated through blood circulation without damaging the surrounding tissues.

THE GLOBUS PROFESSIONAL CAVITATION FREQUENCIES

The cavitation effect is efficiently produced only with low frequencies, lower than 600 KHz, because they are the only one able to penetrate in depth to destroy the adipose tissue. Lower the frequency, deeper is the action of the ultrasonic wave.

- **330 KHz** for Cellulites and medium thickness adiposity (Semi-Professional)
- **526 KHz** to associate the adipose destruction with the production of Noradrenalin (NorAd-526) which helps to burn fat thanks to the increase of the metabolism. Specifically used for the adiposity of the men abdomen
- **3 MHz** for superficial beauty treatments and for skin treatments.

21 PROGRAMS - 30 free programs:

5 ADIPOSE MOBILIZATION NorAd-526
2 CAVITATION 330 KHz
6 BEAUTY
6 SKIN



TECHNICAL FEATURES

Power: 2 W/cm²
 Output: 1 - Head: 1 Ø mm42 – IP65
 Power Supply: main 220V 50/60 Hz or rechargeable batteries (in option)
 Dimensions: 17 x 22 x 6 cm
 Backlight Display
 Emission: continued and pulsed
 Emission Frequency: 3 MHz, 526 KHz e 330 KHz
 CE Device
 Standard: EN 60601-1 EN 60601-1-2

EQUIPMENT

1 carrying bag
 1 ultrasound Lipozero G500
 1 Multi-frequency Head (3 MHz, 526 e 330 KHz)
 1 contact gel bottle
 1 Power pack
 1 operating manual



GLOBUS

LIPOZERO G500

Program List

ADIPOSE MOBILIZATION

NorAd-526

Lipolysis-NorAd-526

Cellulite- NorAd-526

Adiposity- NorAd-526

NorAd-526 1

NorAd-526 2

CAVITATION 330 KHz

Lipolysis

Cellulite

Adiposity abdomen

Orange peel Skin

BEAUTY

Massage

Modelling

Drainage Legs

Drainage arms

Stretch marks

Active principal application

SKIN

Acne

Face Lifting

Wrinkle Treatment

Body spots

Face spots,

Restructuring skin



GLOBUS

GLOBUS ITALIA S.R.L. - Via Vittorio Veneto 52 - 31013 Codognè (TV) - ITALY

Tel (+39)0438.7933 - Fax (+39)0438.793363 - info@globuscorporation.com - www.globuscorporation.com